

**QUALITY FROM THE  
BLACK FOREST**



**IDEAL FOR  
YOUR PRACTICE.**

Dear customers,

are you looking for a medical equipment cart for your practice, dental laboratory or clinic? DEMED is the right choice for you. As the European market leader in medical equipment carts, we manufacture them specifically for use in surgery, endodontics and prophylaxis.

**Still not sure which cart to choose?** On our website [demedmedicalcarts.com](http://demedmedicalcarts.com) you will find all DEMED products. Don't hesitate to contact us for personal, individualized advice.

Have you decided right away? Use the order form on the reverse side. The products you want will be packed securely and delivered ready for use. You can start storing your equipment and tools immediately.

My team and I are happy to be by your side.

With best regards from the Black Forest,



Your sales manager  
Viktor Lebedinski



## MADE IN GERMANY

All DEMED medical carts are developed in Höfen an der Enz in the Black Forest and manufactured in the production facilities of our group of companies.

The carcasses are made of galvanized steel. The handles and standard rails are made of solid, bead-blasted stainless steel. The glass doors are made of ESG safety glass.

You are welcome to arrange a factory tour and convince yourself of the quality and the individual possibilities.

# MEDICAL CARTS THAT SIMPLIFY DAILY LIFE

All DEMED carts are equipped with four retractable ball-bearing casters, doors and drawers with softclose technology.



## E-SERIE

Ample space for accessories, equipment and tools. With a width of only 42 cm, this series takes up very little space.

W x D x H: 420 mm x 500 mm x 900 mm

## P-SERIE

Practical combination of drawers and appliance compartments.

W x D x H: 420 mm x 500 mm x 900 mm

(P-series also available in 800 mm height)

## ADVANTAGES:

- Pure Design
- High-quality components
- Individual configuration
- Wide choice of colors
- Extensive accessories
- Your logo on the glass door



**NEW**

The compact stand-alone cart



3-way socket strip with automatic cable winder integrated in the undershelf

LED-illuminated base, standard colour white

**PX-SERIE**

Our new, super-slim PX model is ideal for housing a single appliance.  
W x D x H: 400 mm x 400 mm x 800 mm

**ERGONOMICS**

Optional standard rail on the drawer to improve ergonomics



**BOTTLE HOLDER**



**DRAWER INSERT**

Made of stainless steel, suitable for various carcasses.



**HOLDER FOR FOOT CONTROL**



Available in three variants



# IN USE AT THE CUSTOMER



„The PX is ideal for small rooms. The visual highlight is the LED-lit base.“



„All devices are ready for immediate use“



„Everything is quickly tidied up again“



Please fill out, photograph and send to [info@demedmedicalcarts.com](mailto:info@demedmedicalcarts.com) or [+49 7081 954550](tel:+497081954550)

Endo Cart E5	Price in €	Colour	Quantity
<b>Endo Cart E5</b>   incl. holder Endomotor (100030) Steel worktop in body colour with all-round stainless steel railing (+3cm)	2.400,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<input type="checkbox"/> <b>Worktop made of solid surface material</b> without railing	50,-	White	
<input type="checkbox"/> <b>Special colour for carcass</b>   Selection from RAL Classic colours	170,-	RAL-Colour: _____	

Endo Cart P1	Price in €	Colour	Quantity
<b>Endo Cart P1</b>   with 4 drawers Steel worktop in body colour with all-round stainless steel railing (+3cm)	1.950,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<b>Endo Cart P1</b>   4 drawers   without pull-out for foot pedal   incl. endo motor holder (100030) incl. standard rail conversion kit (100011) with cable entry (right/left) Steel worktop in body colour with all-round stainless steel railing (+3cm)	2.300,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<b>Endo Cart P1</b>   3 drawers   1 pull-out for foot pedal   incl. endo motor holder (100030) incl. standard rail conversion kit (100011) with cable feed-through (right/left) Steel worktop in body colour with all-round stainless steel railing (+3cm)	2.300,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<input type="checkbox"/> <b>Worktop made of solid surface material</b> without railing	50,-	White	
<input type="checkbox"/> <b>Special colour for body</b>   tion from RAL Classic colours	170,-	RAL-Colour: _____	

Endo Cart PX	Price in €	Colour	Quantity
<b>Endo Cart PX3</b> with 2 drawers   Steel worktop in body colour	1.460,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	

Zubehör	Price in €	Colour	Quantity
<input type="checkbox"/> <b>Stainless steel drawer insert for E5 or P1 incl. four sorting cards</b>	179,-	-	
<input type="checkbox"/> <b>Stainless steel drawer insert for PX</b>	150,-	-	
<input type="checkbox"/> <b>Standard rail at the front of the drawer</b>	70,-	-	
<input type="checkbox"/> <b>Bottle holder</b>	95,-	-	
<input type="checkbox"/> <b>Holder for ultrasound</b>	180,-	-	
<input type="checkbox"/> <b>Holder for foot control with bracket</b>	from 135,-	-	
<input type="checkbox"/> <b>Personalised logo on the glass door for E5 or P1</b>	120,-	-	

All prices are in euros plus VAT and shipping costs. Delivery in cardboard box. The color white corresponds to RAL 9016 and black to RAL 9005 for all medical carts. For further details, please refer to our website demedmedicalcarts.com. The devices shown are not included in the delivery. Subject to technical changes, color variations due to printing and errors.

**Yes, I don't want to miss any more news. Please keep me informed about new DEMED products.**

First name, surname

e-mail

Practice name

Date and signature

Practice stamp

With my signature I agree that my data may be collected and processed for service and quality purposes. As described in the privacy policy at [www.demedmedicalcarts.com](http://www.demedmedicalcarts.com). I have the right to have my personal data updated, changed or deleted at any time with regard to processing and use.

You can find more medical carts and accessories on our website [demedmedicalcarts.com](http://demedmedicalcarts.com)