

**QUALITY FROM THE
BLACK FOREST**



PX3s



PX5s



**IDEAL FOR
YOUR PRACTICE**

SCANNER



Scanner Cart PX3s or PX5s

Ideal for the intraoral scanner in conjunction with a laptop or monitor

- Mobile medical cart made of powder-coated sheet steel
- Dimensions:
 - PX3s: body: 40 cm x 40 cm x 80 cm (W x D x H)
 - PX5s: body: 57 cm x 40 cm x 80 cm (W x D x H)
- 2x drawer with full-extension runner or 1x drawer with full-extension runner and 1x pull-out shelf, multiple socket power strip behind the rear wall with automatic cable winder integrated in the base
- Side-mounted rail system for attaching scanner holder (adapter may be required)
- Brushed stainless steel hanging plate as additional shelf
- White LED strip integrated in the base as standard
- The medical cart is already pre-assembled and is delivered in a cardboard box

List price:

PX3s from 1.560,- EUR

PX5s from 1.590,- EUR

For further questions we are at your disposal. Please contact us for an individual offer.

**Please order by e-mail to: bestellung@demedmedicalcarts.com
or via the order form on the back.**

Your DEMED Sales Team:

Dental Medizintechnik e.K.
Gräfenau 56
D-75339 Höfen an der Enz
bestellung@demedmedicalcarts.com



+49 7081 954550

demedmedicalcarts.com



[DemedMedicalCarts](https://www.instagram.com/DemedMedicalCarts)

Follow us



Please fill out, photograph and send to [✉ bestellung@demedmedicalcarts.com](mailto:bestellung@demedmedicalcarts.com) or [☎ +49 7081 954550](tel:+497081954550)

Medical Cart PX3s	Price in €	Color	Quantity
Mobile cart incl. two H3 drawers (without inserts), side-mounted rail system (+2 cm per side) incl. hanging tray for scanner storage, 3-way socket strip with automatic cable winder and white LED strip integrated in the base	1.560,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<input type="checkbox"/> Special color for body according to RAL classic system	RAL-color: _____		170,-
<input type="checkbox"/> LED lighting RGB (color change) incl. remote control			95,-
<input type="checkbox"/> Stainless steel insert small			150,-
<input type="checkbox"/> Compartment drawer instead of a drawer			w/o surcharge

Medical Cart PX5s	Price in €	Color	Quantity
Mobile cart incl. two H3 drawers, without monitor holder, incl. side-mounted rail system (+2 cm per side), integrated 5-socket power strip with automatic cable winder, white LED strip integrated in the base, door with push-to-open	1.590,-	<input type="checkbox"/> White <input type="checkbox"/> Anthracite	
<input type="checkbox"/> Monitor holder			110,-
<input type="checkbox"/> Special color for body according to RAL classic system	RAL-color: _____		170,-
<input type="checkbox"/> LED lighting RGB (color change) incl. remote control			95,-
<input type="checkbox"/> Stainless steel insert small			150,-
<input type="checkbox"/> Compartment drawer instead of a drawer			w/o surcharge

All prices in Euro plus VAT and shipping costs. Delivery in cardboard box. For all medical carts the color white corresponds to RAL 9016 and anthracite to DB 703. For further details, see our website demedmedicalcarts.com. The devices shown are not included in the delivery. Subject to technical changes, color variations due to printing and errors.

Yes, I don't want to miss any more news. Please keep me informed about new DEMED products.

First name, last name

E-Mail

Name of practice

Date and signature

Practice stamp

With my signature I consent to the recording and processing of my data for service and quality purposes. As listed in the imprint at www.demedmedicalcarts.com. I have the right to obtain at any time the updating, modification or deletion of my personal information regarding its processing and use.

