

**QUALITY FROM THE
BLACK FOREST**

NEWS



Bridge table I5



**IDEAL FOR
YOUR PRACTICE.**



SURGERY

Dear customers,

are you looking for a medical equipment cart for your practice, dental laboratory or clinic? DEMED is the right choice for you. As the European market leader in medical equipment carts, we manufacture them specifically for use in surgery, endodontics and prophylaxis.

Still not sure which cart to choose? On our website demedmedicalcarts.com you will find all DEMED products. Don't hesitate to contact us for personal, individualized advice.

Have you decided right away? Use the order form on the reverse side. The products you want will be packed securely and delivered ready for use. You can start storing your equipment and tools immediately.

My team and I are happy to be by your side.

With best regards from the Black Forest,



Your sales manager
Viktor Lebedinski



MADE IN GERMANY


All DEMED medical carts are developed in Höfen an der Enz in the Black Forest and manufactured in the production facilities of our group of companies.

The carcasses are made of galvanized steel. The handles and standard rails are made of solid, bead-blasted stainless steel. The glass doors are made of ESG safety glass.

You are welcome to arrange a factory tour and convince yourself of the quality and and the individual possibilities.

MEDICAL CARTS THAT SIMPLIFY DAILY LIFE

All DEMED carts are equipped with four retractable ball-bearing casters, doors and drawers with softclose technology.



E-SERIES

Ample space for accessories, equipment and tools. With a width of only 42 cm, this series takes up very little space.

W x D x H: 420 mm x 500 mm x 900 mm

P-SERIES

Practical combination of drawers and equipment compartments.

Available in two widths and two heights.

W x D x H: 420 mm x 500 mm x 900 mm

W x D x H: 420 mm x 500 mm x 800 mm

W x D x H: 500 mm x 500 mm x 800 mm

ADVANTAGES:

- Pure Design
- High-quality components
- Individual configuration
- Wide choice of colors
- Extensive accessories
- Your logo on the glass door



NEWS

Stainless steel
with raised edge
60 x 100 cm

Standard rail
for accessories
included

Stepless electric
height adjustment
from 66 to 125 cm
(without bulky cable)

Stepless height
adjustment from
78 to 120 cm

Shortened base plate
for optimal positioning
on the treatment unit

I-SERIES

The perfect complement to
DEMED carts. High-quality
powder-coated steel frame.
W x D: 650 mm x 450 mm

BRIDGE TABLE

It can be conveniently placed
above the treatment table, and all
instruments are within easy reach.



TABLE TO HOOK

W x D: 400 mm x 480 mm
Can be mounted left or right
Maximum load 5 kg

IN USE BY THE CUSTOMER



"The E-series is also very suitable for small rooms. Thanks to the versatility of the accessories, everything finds its place"



"With the DEMED bridge table, all the tools are always at hand."



"The instrument table is not only practical, but also very elegant"





ORDER

Please fill out, photograph and send to info@demedmedicalcarts.com or [+49 7081 954550](tel:+497081954550)

Medical Carts	Price in €	Color	Quantity
Medical Cart E5 with implant set (100038) Steel worktop with integral stainless steel railing (+3 cm)	2.570,-	<input type="checkbox"/> White	
Medical Cart P1 2 drawers H3, 1 drawer H6 Steel worktop with integral stainless steel railing (+3 cm)	1.880,-	<input type="checkbox"/> White	
<input type="checkbox"/> Worktop in mineral material without railing	50,-	White	
<input type="checkbox"/> Special color for the body Selection from RAL Classic colors	170,-	RAL-Color: _____	

Tables	Price in €	Color	Quantity
Instrument table I3 with standard rail for attachment of accessories	1.250,-	<input type="checkbox"/> White	
Bridge table I5 Stainless steel with raised edge 60 cm x 100 cm, Stepless electric height adjustment, no special color possible	2.050,-	<input type="checkbox"/> White <input type="checkbox"/> Black	
<input type="checkbox"/> Special Color for frame I3 Selection from RAL Classic colors	170,-	RAL-Color: _____	

Accessories	Price in €	Color	Quantity
Table to hook to expand the work surface	390,-	—	
Waste container for standard rail	120,-	—	
Handpiece holder for piezosurgery	179,-	—	
Your logo on the glass door	120,-	—	

All prices are in euros plus VAT and shipping costs. Delivery in cardboard box. The color white corresponds to RAL 9016 and black to RAL 9005 for all medical carts. For further details, please refer to our website demedmedicalcarts.com. The devices shown are not included in the delivery. Subject to technical changes, color variations due to printing and errors.

Yes, I don't want to miss any more news. Please keep me informed about new DEMED products.

First name, last name

E-mail

Name of Practice

Date and signature

Practice stamp

With my signature I consent to the recording and processing of my data for service and quality purposes. As listed in the imprint at www.demedmedicalcarts.com I have the right to obtain at any time the updating, modification or deletion of my personal information regarding its processing and use.

More medical carts and accessories are available on our website demedmedicalcarts.com
We are always available for individual advice.