

**QUALITY FROM THE  
BLACK FOREST**

**NEWS**



Bridge table I5



**IDEAL FOR  
YOUR PRACTICE.**



**SURGERY**

Dear customers,

are you looking for a medical equipment cart for your practice, dental laboratory or clinic? DEMED is the right choice for you. As the European market leader in medical equipment carts, we manufacture them specifically for use in surgery, endodontics and prophylaxis.

**Still not sure which cart to choose?** On our website [demed-medical-carts.de](http://demed-medical-carts.de) you will find all DEMED products-don't hesitate to contact us for personal, individualized advice.

Have you decided right away? Use the order form on the reverse side. The products you want will be packed securely and delivered ready for use. You can start storing your equipment and tools immediately.

My team and I are happy to be by your side.

With best regards from the Black Forest,

Your sales manager  
Viktor Lebedinski

## MAGIC FOR PRACTICE

Scan the QR-code and place the medical cart in your practice.



## MADE IN GERMANY

All DEMED medical carts are developed in Höfen an der Enz, in the Black Forest, and manufactured by ADLER Metal Working - from our group of companies.

The casings are made of galvanized steel. Handles and rails are made of solid and sandblasted stainless steel. Glass doors are made of ESG safety glass.

You are invited to arrange a guided tour of the factory and convince yourself of the quality and possibilities.

# MEDICAL CARTS THAT SIMPLIFY DAILY LIFE

All DEMED carts are equipped with four retractable ball-bearing casters, doors and drawers with softclose technology.



## E-SERIES

Ample space for accessories, equipment and tools. With a width of only 42 cm, this series takes up very little space.

W x D x H: 420 mm x 500 mm x 900 mm

## P-SERIES

Practical combination of drawers and equipment compartments.

Available in two widths and two heights.

W x D x H: 420 mm x 500 mm x 900 mm

W x D x H: 420 mm x 500 mm x 800 mm

W x D x H: 500 mm x 500 mm x 800 mm

## ADVANTAGES:

- Pure Design
- High-quality components
- Individual configuration
- Wide choice of colors
- Extensive accessories
- Your logo on the glass door



**NEWS**



**I-SERIES**

The perfect complement to DEMED carts. High-quality powder-coated steel frame. W x D: 650 mm x 450 mm

**BRIDGE TABLE**

It can be conveniently placed above the treatment table, and all instruments are within easy reach.



**TABLE TO HOOK**

W x D: 400 mm x 480 mm  
Can be mounted left or right  
Maximum load 5 kg

# IN USE BY THE CUSTOMER



"The E-series is also very suitable for small rooms. Thanks to the versatility of the accessories, everything finds its place"



"With the DEMED bridge table, all the tools are always at hand."



"The instrument table is not only practical, but also very elegance"



# ORDER



Please fill out, photograph and send to [info@demed-medical-carts.de](mailto:info@demed-medical-carts.de) or [+49 7081 954550](tel:+497081954550)

Medical Carts	Price in €	Color	Quantity
<b>Medical Cart E5</b>   with implant set (100038) Steel worktop with integral stainless steel railing (+3 cm)	2.449,-	<input type="checkbox"/> White	
<b>Medical Cart P1</b>   2 drawers H3, 1 drawer H6 Steel worktop with integral stainless steel railing (+3 cm)	1.790,-	<input type="checkbox"/> White	
<input type="checkbox"/> <b>Worktop in mineral material</b> without railing	50,-	White	
<input type="checkbox"/> <b>Special color for the carcass</b>   Selection from RAL Classic colors	165,-	RAL-Color: _____	

Tables	Price in €	Color	Quantity
<b>Instrument table I3</b> with standard rail for attachment of accessories	1.250,-	<input type="checkbox"/> White	
<b>Bridge table I5</b> Stainless steel with raised edge 60 cm x 100 cm, Stepless electric height adjustment, <b>no special color possible</b>	1.790,-	<input type="checkbox"/> White <input type="checkbox"/> Black	
<input type="checkbox"/> <b>Special Color for frame I3</b>   Selection from RAL Classic colors	165,-	RAL-Color: _____	

Accessories	Price in €	Color	Quantity
<b>Table to hook</b> to expand the work surface	369,-	—	
<b>Handpiece holder for piezosurgery</b>	170,-	—	
<b>Your logo on the glass door</b>	125,-	—	

All prices are in euros plus VAT and shipping costs. Delivery in cardboard box. White color corresponds to RAL 9016 for all medical carts. For further details, see our website demed-medical-carts.de. The devices shown are not included in the delivery. Subject to technical changes, color variations due to printing and errors.

**Yes, I don't want to miss any more news. Please keep me informed about new DEMED products.**

\_\_\_\_\_  
First name, last name

\_\_\_\_\_  
E-mail

\_\_\_\_\_  
Name of Practice

\_\_\_\_\_  
Date and signature

Practice stamp

With my signature I consent to the recording and processing of my data for service and quality purposes. As listed in the imprint at [www.demed-medical-carts.de](http://www.demed-medical-carts.de). I have the right to obtain at any time the updating, modification or deletion of my personal information regarding its processing and use.

**More medical carts and accessories are available on our website [demed-medical-carts.de](http://demed-medical-carts.de)**  
**We are always available for individual advice.**