

Main catalog 2024



DEMED
MEDICAL CARTS

**MEDICAL CARTS
FOR YOUR PRACTICE**



DemedMedicalCarts

Follow us



+49 7081 954550

demedmedicalcarts.com


MEDICAL CARTS THAT MAKE YOUR DAILY WORK EASIER

Dear DEMED customers,

welcome to the world of DEMED Medical Carts,
your reliable partner for high-quality
medical equipment carts made of steel.

Our medical carts are characterized by their high quality,
functionality and durability. They offer a practical
and efficient solution for the organization and storage of
medical devices, instruments and materials.
In addition, they are ergonomically designed and ease
the everyday work of medical staff.

Discover the variety and quality of DEMED
made in Germany, the leading manufacturer of
medical carts and invest in your medical facility today.



Hüseyin Kanat
Managing director



MADE IN GERMANY

DEMED Medical Carts are developed in Höfen
an der Enz in the Black Forest and
manufactured by ADLER Metal Working -
part of our group of companies.



CONTENTS

DEMED QUALITY **P. 6–13**

Purist design	7
High-quality components	8
Concealed castors & hidden cables	9
Individual configuration & free choice of color	10
How to make your workplace ergonomic	12

PRODUCTS **P. 14–35**

Product overview	15
-------------------------	-----------

Medical carts

E-Series	16
D-Series	18
P-Series	20
PX-Series	22
K-Series	24
S-Series	26

Tables

Instrument table	28
Bridge table	30

Accessories



Which cart suits me?	38
-----------------------------	-----------

SERVICE & ADVICE **P. 42-43**


DEMED QUALITY, CONVINCINGLY GOOD.

Purist, well thought-out designs made from durable materials such as steel, stainless steel or aluminium contribute to the functionality and arrangement of appliances and enable easy cleaning to maintain hygiene standards. Extensive accessories support easy handling in everyday practice.





contribute to the longevity of our medical carts. We only use high-quality components such as blasted handles, standard rails made of stainless steel and doors made of safety glass, this makes our medical carts unique.



CONCEALED CASTORS – AN EMPLOYEE WHO CAN FOLLOW YOU

All DEMED medical carts have smooth-running, integrated swivel castors that enable easy movement and safe transport of your equipment around your premises.

By the way: Our swivel castors come with a lifetime guarantee.



CABLE MANAGEMENT

Integrated cable grommets and a multiple socket with a unique cable retractor are located at key points in the underbody. This makes dangerous tripping hazards and cable breakage a thing of the past.



INDIVIDUAL CONFIGURATION

allows you to customise your medical cart to your own needs and requirements.
Choose from various drawer heights, different drawer pull-outs and storage compartments.
There is space for all devices and instruments.



Fig. shows: E5 series medical carts

FREE CHOICE OF COLOUR

Thanks to a high-quality powder coating in all colours of the RAL Classic coloursystem, your Medical Cart blends in with its surroundings or adds a colour accent if desired.



HOW TO MAKE YOUR WORKPLACE ERGONOMIC

The height adjustability of our instrument and bridge tables means that the work surface can be customised to suit individual requirements. and height of the user adapt.



Integrated castors ensure manoeuvrability and flexibility, which improve work processes and increase efficiency.





Well thought-out accessories are right where you need them.

This promotes an optimum working position and healthy posture.



IMPROVES OUR WORKFLOW

"I am impressed by the size and handling of the bridge table. It has significantly improved treatment procedures and ergonomics. I can adjust the height of the table to my individual needs, which makes my work more comfortable and efficient. I am particularly impressed by the solid workmanship and outstanding durability of the table."

*Prof Angelo Cardarelli
(Master of Piezo Surgery)
from Isernia, Italy*



PRODUCTS




PRODUCT OVERVIEW

MEDICAL CARTS


P. 16-27




E-Series
P. 16



PX-Series
P. 22



D-Series
P. 18



K-Series
P. 24



P-Series
P. 20




S-Series
P. 26

TABLES

P. 28-31




Instrumentent table
P. 28



Bridge table
P. 30


ACCESSORIES


P. 32-37



E-SERIES


Slim cabinet: without taking up valuable space or cluttering up the room, this series impresses with a width of just 42 cm. Yet large enough to store many appliances and slim enough not to take up too much space.





D-SERIES

The series with the largest and highest volume offers plenty of storage space. It is an ideal solution for accommodating wide surgical devices, centrifuges and implant motors.





THE QUALITY IS OUTSTANDING

"As a user from the very beginning, I am still enthusiastic about the many possibilities that the carts and accessory systems from DEMED offer me. It's simply impossible to imagine my practice without the devices and they help me incredibly with efficient endodontic treatment. We haven't had a single problem yet"

Dr. Christoph Sandweg aus Wuppertal



P-SERIES

The professional for prophylaxis, prosthetics and much more. This series offers versatile combination options for a wide range of applications in everyday practice.





THE SPACE PROBLEM WAS SOLVED PERFECTLY

"I'm really impressed with your products! The space problem in my practice is now perfectly solved. The implant motor, piezosurgery and centrifuge now all fit into one cart and are always to hand.

In addition, the ergonomics noticeably improved.

Many thanks for this innovative solution, which has made my working day much easier!"



*Dr Sammy Noubissi
from Washington, USA*

PX-SERIES

Thanks to its compact design with a height of 80 cm the cart takes up very little space and still offers room for laptop, devices and accessories.






LED installed in the underbody



Permanently installed monitor holder on request (PX5s only)

K-SERIES

Orthodontics compact: Your orthodontic accessories, such as pliers and archwires, are stored in our perfectly organized in our Medical Cart with instrument and sorting inserts.





Flexible drawer division




S-SERIES

The cart with the drawers:

Organised storage space - available in two widths, three carcass heights and five drawer heights.

On request with dust protection lips in accordance with RKI guidelines and with optional accessories.





Drawers with full-extension runners and soft-close mechanism



Five drawer heights available

INSTRUMENT TABLE

The mobile instrument table is the perfect addition. With its height-adjustable lift function and the optional standard rail for attaching accessories, it is the ideal assistant.






Optionally with standard rail



BRIDGE TABLE

Wireless, electrically adjustable operating table.
Ensures smooth and jerk-free height adjustment and
can be used for various medical applications.





Button for height adjustment





Stainless steel shelf with raised edge



ACCESSORIES

In addition to our carts, we have numerous solutions such as instrument holders, compartment inserts and much more in our programme.





Holder made of pom for instruments
With or without stainless steel inlay,
autoclavable, suitable for standard rail.



Holder set
For endodontics.



Holder set
For implantology.


Infusion holder
Also available individually.



Pliers holder
Flexible for up to 26 pliers, freely positionable
or for attaching to the standard rail.




Standard rail
Stainless steel solid material, blasted, scope of
delivery incl. spacer and stainless steel screws.



Drawer insert with compartments

Stainless steel insert incl. 4 flexibly positionable
Mouldings, suitable for various carts.



Drawer insert for instruments

Instrument holder incl. 4 strips, max. 52 holders,
Suitable for various carts..



Flexible articulated arm With small standard rail.



Additional tray

For attaching to the standard rail.



Cup holder

For standard cups for
hanging on the standard rail
or stainless steel shelf.




Pull-out tray, tall
Various sizes and colours.



Pull-out tray, flat
Various sizes and colours.



Foot control holder
Stainless steel base plate with foot bracket,
available in different sizes.



Stainless steel tray

To extend the work surface.
Simple attachment to the standard rail.




Contra-angle handpiece holder

Made of stainless steel, for up to 8 contra-angle handpieces.

Litter bin

Frame made of stainless steel.
Simple attachment to the standard rail.



Cylinder lock

For glass doors.



Personalisation

Your logo on the glass door using foil printing.



Dispenser holder

Suitable for up to three dispensing guns, as a table version.




WHICH CART FOR WHICH APPLICATION?

ENDO

SURGERY

STORAGE



E-SERIES

E5 Body (W x D x H): 42 cm x 50 cm x 90 cm

- One drawer with full-extension runner permanently installed
- Three pull-out shelves behind the glass door
- Dimension with standard rail: + 3 cm per side
- Cable outlets: 2 x right, 2 x left, 1 x rear

Worktop:

Optionally steel in body colour with railing or solid surface material in "White"

Equipment:

- Door: Made of frosted toughened safety glass, opens 270°, with steel handle and magnetic catch
- Cable management: 3-way socket strip integrated in the base with automatic cable rewind
- Standard rails: Made of stainless steel, right and Mounted on the left
- Castors: Four smooth-running castors with ball bearings


Colour:

The standard colour is white (Traffic-white RAL 9016). All special colours (RAL Classic) are available at an extra charge.

ENDO

SURGERY

STORAGE



D-SERIES

D1 Carcass (W x D x H): 50 cm x 50 cm x 90 cm

- One drawer with full-extension runner permanently installed
 - Two flat pull-out shelves behind the glass panel
 - Dimension with standard rail: + 3 cm per side
 - Cable outlets: 2 x right, 2 x left, 1 x rear
- D5 Body (W x D x H): 50 cm x 50 cm x 108 cm**
- One drawer with full-extension runner permanently installed
 - Two flat pull-out shelves, one high pull-out shelf
 - Dimension with standard rail: + 3 cm per side
 - Cable outlets: 3 x right, 3 x left, 1 x rear

Worktop:

Optionally steel in body colour with railing or solid surface material in "White"

Equipment:

- Door: Made of frosted toughened safety glass, opens 270°, with steel handle and magnetic catch
- Cable management: 3-way socket strip integrated in the base with automatic cable rewind
- Standard rails: Made of stainless steel, right and Mounted on the left
- Castors: Four smooth-running castors with ball bearings

Colour:


The standard colour is white (Traffic-white RAL 9016). All special colours (RAL Classic) are available at an extra charge.

ENDO

PROPHYLAXE

MOULDING

STORAGE



P-SERIES

P1 carcass (W x D x H): 42 cm x 50 cm x 90 cm

- Three to four drawers with full-extension runners

P3 carcass (W x D x H): 42 cm x 50 cm x 80 cm

- Two to three drawers with full-extension runners

P5 carcass (W x D x H): 50 cm x 50 cm x 80 cm

- Two to three drawers with full-extension runners

Worktop:

Optionally steel in body colour with railing or solid surface material in "White"

Equipment:

- At least one drawer behind the glass door
- Door: Made of frosted toughened safety glass, opens 270°, with steel handle and magnetic catch
- Cable management: 3-way socket strip integrated in the base with automatic cable rewind
- Castors: Four smooth-running castors with ball bearings

Colour:


The standard colour is white (Traffic-whiteRAL 9016). All special colours (RAL Classic) are available at an extra charge.

ENDO

SCANNER

SURGERY

PROPHYLAXE



PX-SERIES

PX3 carcass (W x D x H): 40 cm x 40 cm x 80 cm

Two H3* drawers (without inserts), 3-socket power strip with automatic cable winder and white LED strip integrated in the undershelf

PX3s carcass (W x D x H): 40 cm x 40 cm x 80 cm

Two H3* drawers (without inserts), side-mounted rail system (+2 cm per side) incl. hanging tray for scanner storage, 3-way socket strip with automatic cable rewind and white LED strip integrated in the undershelf

PX5s carcass (W x D x H): 57 cm x 40 cm x 80 cm

Two H3 drawers, without/with monitor holder, incl. side-mounted rail system (+2 cm per side), integrated 5-socket power strip with automatic cable winder, white LED strip integrated in bottom shelf, door with push-to-open function

Equipment:

- Cable outlets: 1 x right, 1 x left, 1 x bottom centre
- Drawers with full-extension runners and shelf pull-outs can be freely combined
- Integrated LED strip: standard colour white (RGB colour change with remote control at an extra charge)
- Castors: Four smooth-running castors with ball bearings

Colour:

The standard colour is white (Traffic-white RAL 9016). All special colours (RAL Classic) are available at an extra charge.

ORTHO- DONTICS



K-SERIES

K1 carcass (W x D x H): 50 cm x 50 cm x 90 cm

- Three drawers with full-extension runners

K3 carcass (W x D x H): 50 cm x 50 cm x 80 cm

- Two drawers with full-extension runners

Worktop:

Choice of steel in body color with railing
or solid surface material in "White"

Equipment:

- Three pull-out shelves behind the glass door
- Door: Made of frosted toughened safety glass, opens 270°, with handle and magnetic catch
- Castors: Four smooth-running castors with ball bearings
- No power supply

Colour:

The standard colour is white
(Traffic-whiteRAL 9016).

All special colours (RAL Classic) are
available at an extra charge.

STORAGE



S-SERIES

S1 carcass (W x D x H): 50 cm x 50 cm x 90 cm

S3 carcass (W x D x H): 42 cm x 50 cm x 90 cm

S5 carcass (W x D x H): 42 cm x 50 cm x 80 cm

S5 B carcass (W x D x H): 50 cm x 50 cm x 80 cm

- Drawers with full-extension runners, in two widths and five heights available
- Castors: Four smooth-running castors with ball bearings
- On request with dust protection lips according to RKI guidelines
- Standard rail optional

Worktop:

Optionally made of steel in body colour with railing or solid surface material in "White"

Colour:

The standard colour is white
(Traffic-whiteRAL 9016).

All special colours (RAL Classic) are
available at an extra charge.


**WE ARE HAPPY TO
ADVISE YOU**

 info@demedmedicalcarts.com



+49 7081 954550

DRIVABLE INSTRUMENT TABLE



I-SERIE

I3 Instrument table

(W x D x H): 65 x 45 cm x 78-120


- Infinitely height-adjustable with lever handle from 78 cm to 120 cm
- Support dimensions: 65 cm x 45 cm
- Material: Frame made of powder-coated steel with polished stainless steel surface
- Castors: Four smooth-running castors with ball bearings

Colour:

The standard colour is white (Traffic-whiteRAL 9016).

All special colours (RAL Classic) are available at an extra charge.

DRIVABLE BRIDGE TABLE



I-SERIE

I5 Bridge table

(W x D x H): 100 x 60 cm x 66-125

- Electrically height-adjustable from 66 cm to 125 cm
- Support dimensions: 60 cm x 100 cm
- Material: Frame made of powder-coated steel with polished stainless steel surface
- Castors: Four smooth-running castors with ball bearings
- Standard rail optional

Colour:

Only available in white (Traffic-whiteRAL 9016) or deep black (RAL 9005)



DemedMedicalCarts

Follow us

www.demedmedicalcarts.com

SERVICE AND ADVICE

You would like to equip your practice, clinic or dental laboratory with our DEMED carts and still need help with the decision-making process, then contact us, we will be happy to advise you.

info@demedmedicalcarts.com

Tel. +49 7081 954550

VISIT OUR
SHOWROOM IN THE
BLACK FOREST



DemedMedicalCarts

Follow us

Notes: The devices and instruments shown in the catalogue are not included in the scope of delivery. Subject to technical changes, printing-related colour deviations and errors. All dimensions, model variants and prices can be found in our separate price list. Photography: PeterOppenländer Fotodesign (www.peter-oppenlaender.de), Other images: Freepik.com, Design: DEMED graphics department



Meet us at a congress or trade fair,
national or international

We will be happy to advise
you by phone, mail,
video call or in your practice



Customised
solutions are adapted
to your needs



All products are
pre-assembled on
delivery and ready
for immediate use

IN-HOUSE PLANNING AND PRODUCTION





DEMED Dental Medizintechnik e.K.
Gräfenau 56 | D-75339 Höfen an der Enz

 info@demedmedicalcarts.com



DemedMedicalCarts

Follow us



+49 7081 954550

demedmedicalcarts.com