

**MADE IN GERMANY**  
QUALITY FROM  
THE BLACK FOREST

**DEMED**  
MEDICAL CARTS



**ARTICLE LIST 2025**

# MEDICAL CARTS **PREMIUM**

INDIVIDUALLY CONFIGURABLE EQUIPMENT  
CARTS OF THE HIGHEST QUALITY FOR ERGONOMIC  
AND EFFICIENT WORK IN YOUR PRACTICE.



# E-SERIES

Body: 42 cm x 50 cm x 90 cm (W x D x H) | dimensions with standard rail: + 3 cm per side

Compartment drawers: 2x flat, 1x high

Cable guides: 2x right, 2x left, 1x back side

Satin glass door, 270° opening, made of toughened safety glass with steel handle and magnetic catch, 3-way socket strip with automatic cable winder and four ball-bearing castors for very smooth running integrated in the base, two stainless steel standard rails (right/left).



Item number	Description/Scope of delivery
100203	<b>Medical Cart E5</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 1x drawer H3*
100202	<b>Medical Cart E5</b> Worktop made of solid surface material in the color „White“, 1x drawer H3*
100210	Compartment drawer (flat) in body color, also available as a retrofit kit
100211	Compartment drawer (high) in body color, also available as a retrofit kit
100060	Stainless steel insert 420 H2* incl. 4 flexibly positionable strips
100061	Stainless steel insert 420 H3* incl. 4 flexibly positionable strips
100077-84	Additional strip for stainless steel drawer inserts (crosswise or lengthwise to match the drawer)
100097	Washable silicone base insert, suitable for drawers and compartment drawers (420)
100009	Special color according to RAL classic system (includes body, compartment drawers and drawers)
100673	Worktop made of solid surface material 420 in special color according to catalog
101119	Extra charge for standard rail on drawer
100045	Personalization: Your logo on the glass door, transfer printing approx 23 cm width, according to customer specifications

\* For more details on drawer heights, general information and shipping costs, see p. 15

# D-SERIES

## Medical Cart D1

Body: 50 cm x 50 cm x 90 cm (W x D x H) | dimensions with standard rail: + 3 cm per side

Compartment drawers: 2x flat

Cable guides: 2x right, 2x left, 1x back side

## Medical Cart D5

Body: 50 cm x 50 cm x 108 cm (W x D x H) | dimensions with standard rail: + 3 cm per side

Compartment drawers: 2x flat, 1x high

Cable guides: 3x right, 3x left, 1x back side

Satin glass door, 270° opening, made of toughened safety glass with steel handle and magnetic catch, 3-way socket strip with automatic cable winder and four ball-bearing castors for very smooth running integrated in the base, two stainless steel standard rails (right/left).



Item number	Description/Scope of delivery
100101	<b>Medical Cart D1</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 1x drawer H3*
100100	<b>Medical Cart D1</b> Worktop made of solid surface material in the color „White“ height difference + 3 cm, 1x drawer H3*
100103	<b>Medical Cart D5</b> Steel worktop in body color with circumferential stainless steel rail, 1x drawer H3*
100102	<b>Medical Cart D5</b> Worktop made of solid surface material in the color „White“, 1x drawer H3*
100110	Additional compartment drawer (flat) in body color, also available as a retrofit kit
100111	Additional compartment drawer (high) in body color, also available as a retrofit kit
100063	Stainless steel insert 500 H2* incl. 4 flexibly positionable strips
100064	Stainless steel insert 500 H3* incl. 4 flexibly positionable strips
100077-84	Additional strip for stainless steel drawer inserts (crosswise or lengthwise to match the drawer)
100098	Washable silicone base insert, suitable for drawers and compartment drawers (500)
100009	Special color according to RAL classic system (includes body compartment drawers and drawers)
100670	Worktop made of solid surface material 500 in special color according to catalog
101119	Extra charge for standard rail on drawer
100045	Personalization: Your logo on the glass door, transfer printing approx 23 cm width, according to customer specifications

\* For more details on drawer heights, general information and shipping costs, see p. 15

# P-SERIES

## Medical Cart P1

Body: 42 cm x 50 cm x 90 cm (W x D x H)

## Medical Cart P3

Body: 42 cm x 50 cm x 80 cm (W x D x H)

## Medical Cart P5

Body: 50 cm x 50 cm x 80 cm (W x D x H)

Cable guides: Rear side top and bottom center

Satin glass door, 270° opening, made of toughened safety glass with steel handle and magnetic catch, 1x flat compartment drawer behind the glass door, height-adjustable, triple socket strip with automatic cable winder and four ball-bearing castors for very smooth running integrated in the base.



Item number	Description/Scope of delivery
100420	<b>Medical Cart P1</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 2x drawer H3*, 1x drawer H6*
100421	<b>Medical Cart P1</b> Worktop made of solid surface material in the color „White“, 2x drawer H3*, 1x drawer H6*
100422	<b>Medical Cart P1</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 4x drawer H3*
100423	<b>Medical Cart P1</b> Worktop made of solid surface material in the color „White“, 4x drawer H3*
100900	<b>Medical Cart P1-Z (Zentrifuge)</b> Steel worktop in body colour with all-round stainless steel railing, height difference + 3 cm, drawers: 1x H4 + 1x H5, 2x pull-out shelves
100901	<b>Medical Cart P1-Z (Zentrifuge)</b> Worktop made of solid surface material, colour 'White', drawers: 1x H4 + 1x H5, 2x pull-out shelves
100450	<b>Medical Cart P3</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 1x drawer H3*, 1x drawer H6*
100452	<b>Medical Cart P3</b> Worktop made of solid surface material in the color „White“, 1x drawer H3*, 1x drawer H6*
100453	<b>Medical Cart P3</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 3x drawer H3*
100454	<b>Medical Cart P3</b> Worktop made of solid surface material in the color „White“, 3x drawer H3*
100455	<b>Medical Cart P5</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 3x drawer H3*
100456	<b>Medical Cart P5</b> Worktop made of solid surface material in the color „White“, 3x drawer H3*
100011	Standard rail conversion kit: 1x cable guide + 1x stainless steel standard rail mounted right or left
100009	Special color according to RAL classic system (includes body, compartment drawers and drawers)
100673	Worktop made of solid surface material 420 in special color according to catalog
100670	Worktop made of solid surface material 500 in special color according to catalog
101119	Extra charge for standard rail on drawer

# PX-SERIES

**Medical Cart PX3**

Body: 40 cm x 40 cm x 80 cm (W x D x H)

**Medical Cart PX3s**

Body: 40 cm x 40 cm x 80 cm (W x D x H)

**Medical Cart PX5s**

Body: 57 cm x 40 cm x 80 cm (W x D x H)

Cable guides: 1x right, 1x left, 1x bottom center

Drawers with full-extension runners and compartment drawers can be freely combined, socket strip with automatic cable winder and four ball-bearing castors for very smooth running integrated in the base.



Item number	Description/Scope of delivery
101123	<b>Medical Cart PX3</b> Mobile cart incl. two H3* drawers (without inserts), 3-socket power strip with automatic cable winder and white LED strip integrated in the base
101124	<b>Medical Cart PX3s</b> Mobile cart incl. two H3* drawers (without inserts), side-mounted rail system (+2 cm per side) incl. hanging tray for scanner storage, 3-way socket power strip with automatic cable winder and white LED strip integrated in the base
101144	<b>Medical Cart PX5s without monitor holder</b> Mobile cart incl. two H3* drawers, without monitor holder, incl. side-mounted rail system (+2 cm per side), integrated 5-socket power strip with automatic cable winder, white LED strip integrated in the base, door with push-to-open
101145	<b>Medical Cart PX5s with monitor holder</b> Mobile cart incl. two H3* drawers, with monitor holder, incl. side-mounted rail system (+2 cm per side), integrated 5-socket power strip with automatic cable winder, white LED strip integrated in the base, door with push-to-open
100009	Special color according to RAL classic system (includes body, compartment drawers and drawers)
100067	Stainless steel insert small
101146/-47	Compartment drawer instead of a drawer

\* General information and shipping costs, see p. 15

# K-SERIES

## Medical Cart K1

Body: 50 cm x 50 cm x 90 cm (W x D x H)

## Medical Cart K3

Body: 50 cm x 50 cm x 80 cm (W x D x H)

Compartment drawer: 1x flat, 2x high, 1x instrument holder,

1x flexible drawers plier holder, 2x stainless steel inserts with 6 flexibly positionable strips.

Satin glass door, 270° opening, made of toughened safety glass with steel handle and magnetic catch, four ball-bearing castors for very smooth running integrated in the base.



Item number	Description/Scope of delivery
100822	<b>Medical Cart K1</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 1x drawer H5*, 1x drawer H3*, 1x drawer H4*
100823	<b>Medical Cart K1</b> Worktop made of solid surface material in the color „White“, 1x drawer H5*, 1x drawer H3*, 1x drawer H4*
100824	<b>Medical Cart K3</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm, 1x drawer H4*, 1x drawer H5*
100825	<b>Medical Cart K3</b> Worktop made of solid surface material in the color „White“, 1x drawer H4*, 1x drawer H5*
100013	Plier holder of stainless steel
100056	Instrument holder 500 incl. 4 flexibly positionable strips, max. 52 holders
100009	Special color according to RAL classic system (includes body, compartment drawers and drawers)
100670	Worktop made of solid surface material 500 in special color according to catalog
100098	Bottom insert made of silicon 500, washable, suitable for compartment and drawers
100045	Personalization: Your logo on the glass door, transfer printing approx 23 cm width, according to customer specifications

\* For more details on drawer heights, general information and shipping costs, see p. 15

# S-SERIES

## Medical Cart S1

Body: 50 cm x 50 cm x 90 cm (W x D x H)

## Medical Cart S3

Body: 42 cm x 50 cm x 90 cm (W x D x H)

## Medical Cart S5

Body: 42 cm x 50 cm x 80 cm (W x D x H)

## Medical Cart S5 B

Body: 50 cm x 50 cm x 80 cm (W x D x H)



Body without drawers, four ball-bearing castors for very smooth running integrated in the base.

### Choose your individual drawer combination from five drawer heights. (See below)

Item number	Description/Scope of delivery
100401	<b>Medical Cart S1 height units (H): 24</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm
100400	<b>Medical Cart S1 height units (H): 24</b> Worktop made of solid surface material in the color „White“
100403	<b>Medical Cart S3 height units (H): 24</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm
100404	<b>Medical Cart S3 height units (H): 24</b> Worktop made of solid surface material in the color „White“
100405	<b>Medical Cart S5 height units (H): 21</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm
100406	<b>Medical Cart S5 height units (H): 21</b> Worktop made of solid surface material in the color „White“
101015	<b>Medical Cart S B height units (H): 21</b> Steel worktop in body color with circumferential stainless steel rail, height difference + 3 cm
101016	<b>Medical Cart S B height units (H): 21</b> Worktop made of solid surface material in the color „White“
100407	Dust sealing lips according to the RKI: All drawer fronts with sealing lip and individual bars in body color attached to the body between the drawers
100060/-63	Stainless steel insert 420/500 H2* incl. 4 flexibly positionable strips
100061/-64	Stainless steel insert 420/500 H3* incl. 4 flexibly positionable strips
100077-84	Additional strip for stainless steel drawer inserts (crosswise or lengthwise to match the drawer)
100097/	
100098	Washable silicone base insert, suitable for drawers and compartment drawers (420/500)
100009	Special color according to RAL classic system (includes body, compartment drawers and drawers)

### Drawers

Choose from five drawer heights from H2 to H6, also in combination! A total of 24 height units (H) are available for S1 and S3, and 21 height units (H) for S5 and S5B. All drawers of the S-series are equipped with a full-extension runner.

e.g. 4x H3 + 2x H4 + 2x H2 = 24 height units

Height units	H2	H3	H4	H5	H6
Internal height (mm)	34	66	98	130	162
Height of cover (mm)	61	93	125	157	189



# TRAY-CART

Body: 50 cm x 50 cm x 90 cm (W x D x H)

four ball-bearing castors for very smooth running  
integrated in the base


Space for up to 22 standard trays

Flexible stainless steel holder system

Medium space partitioning (dust-tight)

Steel worktop in body colour

Aluminium doors in body colour



Item number	Description/Scope of delivery
100410	<b>Tray-Cart</b> incl. 44 stainless steel holders (for up to 22 standard trays)
100009	Special color according to RAL classic system

General information and shipping costs see p. 15


# I-SERIES

## Instrument table I3

Infinitely height-adjustable from 78 cm to 120 cm

Support dimensions: 65 cm x 45 cm

Four ball-bearing castors for very smooth running integrated in the base, frame made of powder-coated steel, polished stainless steel support with raised edge, shortened base plate for optimum positioning on the treatment unit.




## Bridge table I5

Electrically height-adjustable from 66 cm to 125 cm

Support dimensions: 60 cm x 100 cm

Four ball-bearing castors for very smooth running integrated in the base, frame made of powder-coated steel, polished stainless steel support with raised edge available in RAL 9016 traffic white or RAL 9005 black.



Item number	Description/Scope of delivery
100352	<b>Instrument table I3</b> without standard rail
100353	<b>Instrument table I3</b> incl. standard rail for attaching accessories
100354	<b>Bridge table I5</b> Mobile bridge table made of powder-coated steel, shelf made of polished stainless steel with raised edge 60 cm x 100 cm, wireless electric height adjustment (Li-ion battery) from 66 cm - 125 cm by means of push-button
100053	Cup holder, for standard cups for attaching on standard rail or small stainless steel tray
100017	Waste bin, frame stainless steel incl. hanging device for standard rail (without bag)
100009	Special color according to RAL classic system (only applies to I3)

General information and shipping costs see p. 15

# ACCESSORIES **PREMIUM**

**INSTRUMENT HOLDERS, COMPARTMENT  
INSERTS AND MUCH MORE.**



# DEMED ACCESSORIES

## HOLDER FOR STRAIGHT AND CONTRA-ANGLE HANDPIECES

Our holders are made of POM (polyoxymethylene, thermoplastic material) in the standard sizes of the hand-held and hand-operated versions and contra-angle handpieces. These are suitable for the standard rails we fit and can be cleaned in the practice without hesitation using a thermal disinfectant. The stainless steel inlays are bead-blasted in accordance with MPG stainless steel. This gives them a durable and matt texture.

Item number	Description/Scope of delivery
100030	Holder with stainless steel inlay Ø 22/26, suitable for implant and endomotor, autoclavable, suitable for standard rail
100031	Holder with stainless steel inlay Ø 18/22, suitable e.g. for piezotomes, autoclavable, suitable for standard rail
10003-x	Straight and contra-angle handpiece holder, without stainless steel inlay, available in sizes Ø 12, 15, 18, 22 mm
100037	Holder set Endo, 1x handpiece holder with stainless steel inlay Ø 22/26 mm, 1x holder POM Ø 18 mm, 1x holder POM Ø 22 mm
100038	Implant holder set, 1x handpiece holder with stainless steel inlay Ø 22/26 mm, 1x Infusion holder POM incl. stainless steel gallows
100039	Holder for laser, two mounts in one holder: Ø 10 mm + Ø 12 mm
100042	Holder as special size, manufactured according to customer specifications (max. size according to standard holder 30 mm)
100020	stainless steel boom (holds infusion bottle)
100021	POM holder for stainless steel boom
100018	Holder for dispenser, suitable for up to three dispensing guns
101019	8-position contra-angle handpiece holder made of matt and blasted stainless steel incl. holder made of POM

General information and shipping costs see p. 15

## DRAWER INSERTS AND EXTENSIONS

Item number	Description/Scope of delivery
100060	Stainless steel insert 420 H2* incl. 4 flexibly positionable strips (suitable for cabinet width 42 cm)
100061	Stainless steel insert 420 H3* incl. 4 flexibly positionable strips (suitable for cabinet width 42 cm)
100063	Stainless steel insert 500 H2* incl. 4 flexibly positionable strips (suitable for cabinet width 50 cm)
100064	Stainless steel insert 500 H3* incl. 4 flexibly positionable strips (suitable for cabinet width 50 cm)
100077-84	Additional strip for stainless steel drawer inserts (crosswise or lengthwise to match the drawer)
100055	Instrument holder 420 incl. 4 flexibly positionable strips, max. 52 holders (suitable for cabinet width 42 cm)
100056	Instrument holder 500 incl. 4 flexibly positionable strips, max. 52 holders (suitable for cabinet width 50 cm)
100110	Additional compartment drawer 500, flat
100111	Additional compartment drawer 500, high
100210	Additional compartment drawer 420, flat
100211	Additional compartment drawer 420, high
100673	Worktop made of solid surface material 420 in special color according to catalog
100670	Worktop made of solid surface material 500 in special color according to catalog
100097	Washable silicone base insert, suitable for drawers and compartment drawers (420)
100098	Washable silicone base insert, suitable for drawers and compartment drawers (500)

General information and shipping costs see p. 15

# STANDARD RAIL AND EXTENSIONS

Item number	Description/Scope of delivery
100013	Plier holder, flexible for up to 26 pliers, freely positionable or for attaching to the standard rail
100014	Holder for foot control, stainless steel base plate with foot bracket, small (for one foot control)
100019	Holder for foot control, stainless steel base plate with foot bracket, round (for a foot control, e.g. Piezotome)
100015	Holder for foot control, stainless steel base plate with foot bracket, large (for up to 3 foot controls)
100052	Stainless steel shelf for accessories and appliances, for simple and practical hanging on the standard rail
100053	Cup holder for standard cups for attaching on standard rail or small stainless steel tray
100017	Waste bin, frame stainless steel incl. hanging device for standard rail (without bag)
100010	Standard rail, stainless steel solid material, blasted, scope of delivery incl. spacer and stainless steel screws M6
100016	Stainless steel tray to extend the work surface, easy to attach to the standard rail
101020	Flexible articulated arm with small standard rail for up to two holders

General information and shipping costs see p. 15

# GENERAL INFORMATION ON ALL MEDICAL CARTS PREMIUM AND TABLES

**Color:** All our medical carts and instrument tables are powder-coated in the color RAL 9016 traffic white, unless otherwise ordered.

**Special colors:** We can produce all colors of the RAL-Classic color system in our powder coating.

The surcharge for special colors is € 170.00 and includes all components. If several medical carts or instrument tables are ordered in one order, the surcharge is only charged once. For color combinations (e.g. different colors of body and drawer) we create an individual offer.

**Power supply:** Our D-, E-, PX- and P-series medical carts are equipped with an integrated power supply. This includes a socket strip on the rear wall on the inside. D, E, PX3 and P series: 3-way socket strip, PX5: 5-way socket strip I5: wireless height-adjustable (Li-ion battery) via push-button (alternatively height-adjustable via power cable)

The connection is made using an extendable extension cable with a length of approx. three meters.

**Doors:** The doors are made of satinized safety glass (ESG) as standard. Opening to the left as standard. Can also open to the right on request.

## SHIPPING COSTS

Description/Scope of delivery

Medical carts pre-assembled, bridge table

each Instrument table & PX series pre-assembled

Parts and accessories, up to 3 kg

Parts and accessories, from 3 kg

**Shipping costs:** Shipping costs are requested from the forwarding agents and invoiced to the customer.

**Note:** Delivery by forwarding agent. Goods are not carried into the practice.


**General:** All prices are in euros, plus statutory VAT. All our medical carts are pre-assembled. Illustrations and information on brochures and price lists are for guidance only. Errors and omissions excepted. The general terms and conditions of DEMED Dental Medizintechnik e.K. apply, see p. 22-23.

## DRAWERS

Choose from five drawer heights from H2 to H6, also in combination! A total of 24 height units (H) are available for S1 and S3, or 21 height units (H) for S5 and S5B. All drawers in the S-series are equipped with a full-extension runner.

e.g. 1x H2 + 2x H4 + 2x H3 + 1x H5 = 21 height units

Height units	H2	H3	H4	H5	H6
Internal height (mm)	34	66	98	130	162
Height of cover (mm)	61	93	125	157	189



# MEDICAL CARTS **BASIC**

HIGH-QUALITY MEDICAL MULTI-PURPOSE CARTS  
IDEAL AS A SIDE TABLE OR AS AN EQUIPMENT  
STORAGE UNIT IN YOUR PRACTICE.





# DL-SERIES

## DL24-40

Body: 40 cm x 45 cm x 73 cm (W x D x H)

Drawer inside: 30,2 cm x 38,4 cm x 74 cm (W x D x H)

## DL24-60

Body: 60 cm x 45 cm x 73 cm (W x D x H)

Drawer inside: 50,2 cm x 38,4 cm x 74 cm (W x D x H)

Frame with all-round railing for a secure hold

Drawer with full extension and soft close technology


4x swivel castors, two with brakes

Frame made of sheet steel with high-quality powder coating

2 cable apertures for easy cable routing

Max. Load capacity: 50 kg; Max. Max. load capacity drawer: 10 kg

Simple assembly for quick use



Item number	Description/Scope of delivery
100380	<b>DL24-40</b> White colour
100381	<b>DL24-40</b> Black colour
100381	<b>DL24-40</b> Special color according to RAL classic system
100383	<b>DL24-60</b> White colour
100384	<b>DL24-60</b> Black colour
100385	<b>DL24-60</b> Special color according to RAL classic system

General information and shipping costs see p. 20

# CX-SERIES

Body: 40 cm x 38 cm x 85 cm (W x D x H)

4x swivel castors, two with brakes  
 Upper shelf with all-round railing for a secure hold  
 2 cable outlets for neat cable routing  
 Base plate with 3 openings for foot pedals and air and water accessories  
 Integrated triple socket with mains cable  
 Material: Sheet steel, powder-coated in white  
 Simple assembly for quick use



Item number	Description/Scope of delivery
100375	<b>CX</b> White colour
100376	<b>CX</b> Black colour
100377	<b>CX</b> Special color according to RAL classic system

General information and shipping costs see p. 20

# S24-SERIES

Body: 57 cm x 43 cm x 85 cm (W x D x H)

4x swivel castors, two with brakes

Integrated cable feed-through in the table top

3-way socket with 3 m supply cable, integrated in the frame to save space


Cart can be loaded up to 10 kg

Optionally available on request:

Standard rail mountable on both sides

Holder for monitor, power supply unit or power bank can be retrofitted

Personalisation with your own logo



Item number	Description/Scope of delivery
101149	<b>S24</b> White colour
101150	<b>S24</b> Black colour
101152	<b>S24</b> Special color according to RAL classic system

General information and shipping costs see p. 20

# GENERAL INFORMATION ON ALL MEDICAL CARTS **BASIC**

**Colour:** Unless otherwise ordered, all our medical carts are powder-coated in the colour RAL 9016 traffic white. powder-coated. The colour black corresponds to RAL 9005.

**Special colours:** We can produce all colours of the RAL Classic colour system in our powder coating. The surcharge for special colours is € 85.00. If you order several medical carts or instrument tables in one order, the surcharge is only charged once.

**Power supply:** The CX series is equipped with an integrated power supply. This includes a 3-socket power strip on the rear wall on the inside. The connection is made using a three metre power cable. The sockets comply with the Euro standard (Schuko). For deliveries to Switzerland or England, these are supplied in accordance with the country code (main cable by means of an adapter, which is included).

---

## SHIPPING COSTS

Description/Scope of delivery

---

Medical Carts Basic, pre-assembled elements, supplied in a box

---

**Shipping costs:** Shipping costs apply nationwide (except islands). For larger quantities or deliveries abroad, the price is on request.

**General:** All prices are in euros, plus statutory VAT. Medical Carts Basic are delivered as pre-assembled elements, packed in a cardboard box. Delivery is made by delivery service Invoice amounts are due within 7 days of the invoice date without deduction. Illustrations and information on brochures and price lists are for guidance only. We reserve the right to make changes and errors. The General Terms and Conditions of General Terms and Conditions of DEMED Dental Medizintechnik e.K. apply, see p. 22-23.

# NOTES FOR MY ORDER:

---

---

---

---

---

Item number      Description/Scope of delivery

---

---

---

---

---

---

---

---

---

---

---



 [bestellung@demedmedicalcarts.com](mailto:bestellung@demedmedicalcarts.com)

 [demedmedicalcarts.com](http://demedmedicalcarts.com)

 **+49 7081 954550**

# GENERAL TERMS AND CONDITIONS

## 1. Validity of the conditions

- a) The deliveries, services and offers of the seller are made exclusively on the basis of these terms and conditions. They shall therefore also apply to all future business relations, even if they are not expressly agreed again. These terms and conditions shall be deemed accepted at the latest upon receipt of the goods or services. Counter-confirmations by the buyer with reference to his terms and conditions of business or purchase are hereby rejected.
- b) All agreements made between the seller and the buyer for the purpose of executing this contract must be recorded in writing.

## 2. Offer and conclusion of contract

- a) The Seller's offers are subject to change and non-binding. Declarations of acceptance and all orders require the written confirmation of the seller to be legally effective.
- b) Drawings, illustrations, dimensions, weights or other performance data are only binding if this has been expressly agreed in writing.
- c) The seller's sales employees/agents are not authorized to make verbal collateral agreements or give verbal assurances that go beyond the content of the written contract.
- d) Custom-made products and special colors cannot be exchanged, except in the case of defects.

## 3. Prices

- a) Unless otherwise stated, the Seller shall be bound by the prices contained in its offers for 14 days from their date. Otherwise, the prices stated in the Seller's order confirmation plus the respective statutory value added tax shall apply. Additional deliveries and services shall be invoiced separately.
- b) Unless otherwise agreed, prices are quoted in euros plus VAT, packaging, freight and postage.

## 4. Delivery and performance time

- a) Delivery dates or deadlines that have been agreed as binding or non-binding are can be made must be in writing.
- b) Delays in delivery and performance due to force majeure and due to events that make it impossible for the seller to make delivery.

The Seller shall not be responsible for any circumstances that temporarily make delivery or performance significantly more difficult or impossible - including, in particular, strikes, lockouts, official orders, etc., even if they occur at the Seller's suppliers or their subcontractors - even in the case of bindingly agreed deadlines and dates. They entitle the seller to postpone the delivery or service for the duration of the hindrance plus a reasonable start-up time or to withdraw from the contract in whole or in part due to the part not yet fulfilled.

- c) If the hindrance lasts longer than three months, the buyer is entitled, after setting a reasonable grace period, to withdraw from the contract with regard to the part not yet fulfilled. If the delivery time is extended or if the seller is released from his obligation, the buyer cannot derive any claims for damages from this. The Seller may only invoke the aforementioned circumstances if he notifies the Buyer immediately.
- d) If the seller can foresee that the goods cannot be delivered on time, the business partner will be informed immediately and in writing, the reasons for this will be communicated to him and, if possible, the expected delivery date will be stated.

- e) The Seller shall be entitled to make partial deliveries and render partial services at any time, unless the partial delivery or partial service is of no interest to the Buyer.
- f) Compliance with the Seller's delivery and performance obligations requires the timely and proper fulfillment of the Buyer's obligations.
- g) If the Buyer is in default of acceptance, the Seller shall be entitled to demand compensation for the damage incurred; the risk of accidental deterioration and accidental loss shall pass to the Buyer upon the occurrence of default of acceptance.

## 5. Transfer of risk

The risk shall pass to the Buyer as soon as the consignment has been handed over to the person carrying out the transportation or has left the Seller's warehouse for the purpose of shipment. If shipment is delayed at the buyer's request, the risk shall pass to the buyer upon notification of readiness for shipment.

## 6. Rights of the buyer due to defects

- a) The products are delivered free of manufacturing and material defects; the period for asserting claims for defects is one year from delivery of the products.
- b) If the seller's operating or maintenance instructions are not followed, changes are made to the products, parts are replaced or consumables are used that do not comply with the original specifications, claims for defects in the products shall lapse if the buyer does not refute a substantiated claim that one of these circumstances caused the defect.
- c) The Buyer must notify the Seller's customer service department of defects in writing without delay, but at the latest within one week of receipt of the delivery item. Defects that cannot be discovered even after careful inspection within this period must be reported to the seller in writing immediately after discovery.
- d) In the event that the Buyer notifies the Seller that the Products are defective, the Seller shall, at its option and expense, require that:
  - the defective part or device is sent to the seller for repair and subsequent return;
  - the buyer keeps the defective part or device ready and a service technician of the seller is sent to the buyer to carry out the repair.

If the Buyer requests that rectification work be carried out at a location specified by him, the Seller may comply with this request, whereby replaced parts shall not be charged, while working time and travel expenses shall be paid at the Seller's standard rates.

- e) If the rectification of defects fails after a reasonable period of time, the buyer may, at his discretion, demand a reduction in payment or withdraw from the contract.
- f) Liability for normal wear and tear is excluded.
- g) Claims for defects against the seller are only available to the direct purchaser and are not assignable.

## 7. Retention of title

- a) Until the fulfillment of all claims (including all current account balances) to which the seller is entitled against the buyer now or in the future for any legal reason, the seller is granted the following securities, which he will release on request at his discretion, insofar as their value exceeds the claims by more than 20%.
- b) The goods remain the property of the seller. Processing or transformation shall always be carried out for the seller as manufacturer, but without any obligation for him. If the seller's (co-)ownership expires due to combination, it is hereby agreed that the buyer's (co-)ownership of the uniform item shall pass to the seller in proportion to its value (invoice value). The Buyer shall store the Seller's (joint) property free of charge. Goods to which the seller is entitled to (co-)ownership are hereinafter referred to as reserved goods.
- c) The buyer is entitled to process and sell the reserved goods in the ordinary course of business as long as he is not in default. Pledges or transfers by way of security are not permitted. The Buyer hereby assigns to the Seller in full by way of security any claims arising from the resale or any other legal grounds (insurance, tort) in respect of the goods subject to retention of title (including balance claims arising from a current account). The seller revocably authorizes him to collect the claims assigned to the seller for the seller's account in his own name. This collection authorization can only be revoked if the buyer does not properly meet his payment obligations.
- d) In the event of access by third parties to the reserved goods, in particular seizures, the Buyer shall draw attention to the Seller's ownership and inform the Seller immediately so that the Seller can enforce its ownership rights. Insofar as the third party is not in a position to reimburse the seller for the costs incurred in this connection, the seller shall judicial or extra-judicial costs, the person liable for these costs is the buyer.
- e) In the event of breach of contract by the Buyer - in particular default in payment - the Seller shall be entitled to withdraw from the contract and demand the return of the reserved goods.

## 8. Payment

- a) Payment is made in accordance with the customer agreement. The Seller shall be entitled, despite the Buyer's provisions to the contrary, to offset payments against the Buyer's older debts first and shall inform the Buyer of the type of offsetting that has taken place. If costs and interest have already been incurred, the seller shall be entitled to offset the payment first against the costs, then against the interest and finally against the principal performance.
- b) Payment shall only be deemed to have been made when the seller can dispose of the amount. In the case of checks, payment shall not be deemed to have been made until the check has been cashed.
- c) If the Buyer is in default, the Seller shall be entitled to demand interest at a rate of 12.5% above the base rate as liquidated damages from the relevant date. They are to be set lower if the buyer proves a lower charge; proof of higher damages by the seller is permissible.
- d) If the seller becomes aware of circumstances that call into question the creditworthiness of the buyer, in particular if a check is not honored or payments are suspended, or if the seller becomes aware of other circumstances that call into question the creditworthiness of the buyer, the seller is entitled to declare the entire remaining debt due, even if he has accepted checks. In this case, the seller is also entitled to demand advance payments or security deposits.

- e) The buyer is only entitled to offset, withhold or reduce payment, even if complaints or counterclaims are asserted, if the counterclaims have been legally established or are undisputed. However, the buyer is also entitled to withhold payment due to counterclaims arising from the same contractual relationship.

## 9. Design changes

The seller reserves the right to make design changes at any time; however, he is not obliged to make such changes to products that have already been delivered.

## 10. Liability

- a) Claims for damages are excluded irrespective of the type of breach of duty, including tortious acts, except in cases of willful misconduct or gross negligence.
- b) In the event of a breach of material contractual obligations, the Seller shall be liable for any negligence, but only up to the amount of the foreseeable damage. Claims for loss of profit, saved expenses, claims for damages from third parties and other indirect and consequential damages cannot be demanded unless a quality feature guaranteed by the seller is specifically intended to protect the buyer against such damages.
- c) The limitations and exclusions of liability in the paragraphs a) and b) do not apply to claims that have arisen due to fraudulent behavior on the part of the seller, nor to liability for guaranteed characteristics, for claims under the Product Liability Act or for damages resulting from injury to life, limb or health.
- d) Insofar as the seller's liability is excluded or limited, this shall also apply to employees, workers, representatives and vicarious agents of the seller.

## 11. Applicable law, place of jurisdiction, partial invalidity

- a) The law of the Federal Republic of Germany shall apply to these terms and conditions and the entire legal relationship between the seller and the buyer.
- b) If the purchaser is a merchant or a legal entity under public law, Calw shall be the exclusive place of jurisdiction for all disputes arising directly or indirectly from the contractual relationship.
- c) Should a provision in these terms and conditions or a provision within the framework of other agreements be or become invalid, this shall not affect the validity of all other provisions or agreements.

Höfen an der Enz, February 2025



DEMED Dental Medizintechnik e.K.  
Gräfenau 56 | D-75339 Höfen an der Enz

 [info@demedmedicalcarts.com](mailto:info@demedmedicalcarts.com)

 [demedmedicalcarts.com](http://demedmedicalcarts.com)

 [+49 7081 954550](tel:+497081954550)

   
DemedMedicalCarts  
Follow us