




TRAY-CART

Mobile equipment trolley made of steel

- Dimensions: 50 cm x 50 cm x 90 cm (W x D x H)
- Space for up to 22 standard trays
- Flexible stainless steel holder system
- Medium space partitioning (dust-tight)
- Four smooth-running castors with ball bearings
- Steel worktop in body colour
- Standard colour white (RAL 9016)
- Special colour (RAL-Classic) available for an extra charge (175,- Euro)

List price: 1.550,- Euro

Art.-Nr. 100410



We are happy to answer any further questions you may have. Please contact us for a customised offer. Please send your order by e-mail to: bestellung@demedmedicalcarts.com

Your DEMED Sales Team:

Dental Medizintechnik e.K.

Gräfenau 56

D-75339 Höfen an der Enz

bestellung@demedmedicalcarts.com



+49 7081 954550

demedmedicalcarts.com



[DemedMedicalCarts](https://www.instagram.com/DemedMedicalCarts)

Follow us